



Fonctions de base:

	Back up		Relax position
	Back down		Sleep position
	Knee up		Nurse position <i>only HD</i>
	Knee down		Anti-Trendelenburg
	HL up		UBL
	HL down		Trendelenburg <i>only HD</i>
			Lock function <i>only HD</i>

Fonctions supplémentaires:

	-L-		Lock All		SAVE Transport position
	-R-		Unlock All (when error = Reset)		Trendelenburg HD handset (magnetic lock)
	-C-		CPR		
	-T-		Transport position		